

NT Surgical Kit

Note: * means optional

059.0460



Precision drill



059.0020T

Twist drill

L=7-15mm



Φ=2.0mm

059.0030T

Twist step drill

L=7-15mm



Φ2.4/2.8mm

059.0060T

Φ2.8/3.2mm

059.0090T

Φ3.2/3.6mm

059.0120T

Φ3.8/4.2mm

059.0150T

Φ4.2/4.6mm

059.0180T

Drill extension



026.0086

Parallel pin



Φ=2.0mm

059.0300

Connector



059.0410

Tap

L=7-15mm



Φ 3.5mm

* 059.0230T

Φ 4.3mm

* 059.0240T

Φ 5.0mm

* 059.0250T

Screwdriver



059.0480



Handpiece

L=27mm

L=32mm

059.0440

059.0450

Implant driver

L=28mm

L=37mm



NP

059.0370

059.0371



RP

059.0380

059.0381

Multi-unit abutment driver



L=19mm

038.1383

Torque control device



034.0570

Depth gauge



059.0310

Screwdriver NT



L=21mm

* 059.0580